

SIMPLICITY WORKFLOW

Beautifully Simple, from Prep to Polish



Kerr[™]

INDIRECT RESTORATION

STEP 1 - PREP



Indirect Preparation

STEP 2 - IMPRESSION



Take 1™ Advanced™

STEP 4 - FINAL RESTORATION



Maxcem Elite™ Universal & Maxcem Elite™ Chroma

STEP 3 - TEMPORARY SOLUTION



TempSpan™ & Temp-Bond™

STEP 5 - POLISH



Identiflex™ Diamond Ceramic Polishers

Twin-zone diamond polishers for treating ceramics chairside.



OptiShine™

Concave-shaped polishing brush with embedded polishing particles for composite, compomer, resin-modified glass ionomer or ceramics.

DIRECT RESTORATION

STEP 1 - PREP



Cavity Preparation

STEP 2 - BOND



OptiBond™ Universal

STEP 4 - CURE



Demi™ Plus Curing Light

STEP 3 - RESTORE



SimpliShade™ Universal



SimpliShade™ Bulk Fill



SimpliShade™ Bulk Fill Flow

STEP 5 - POLISH



OptiDisc™

Transparent and flexible discs in 4 different grits and 3 diameters for composites, glass ionomers, amalgams, semi and precious metals.



Opti1Step™

Multiple use composite polishing system in 3 shapes with embedded polishing particles.



Occlubrush™

One-step polishing system. Special fibers with in-build Silicone-carbide abrasive particles. No paste is necessary.



OptiBond™ Universal

One Adhesive, Many Procedures

-  **Exceptional Bonding** - Built-in GPDM monomer provides excellent bond strength to all surfaces and substrates for indirect restorations.
-  **Simplified Steps** - Streamlined One Bottle System and one-coat application eliminate steps for all bonding techniques including self-etch.



SimpliShade™

Simplified Universal Composite System

-  **Shade-Matching Made Easy** - With just three shades (Light, Medium, Dark) to match all 16 VITA® classical shades, matching much faster and easier.
-  **Universal For All Cavity Classes** - Use SimpliShade with confidence in all cavity classes without the need for a blocker



Maxcem Elite™ Universal and Maxcem Elite™ Chroma

One Cement for All Your Needs

-  **Freedom to Choose** - Simplify your inventory requirements and keep using your preferred cementation technique with a self-adhesive cement that can also be used as a resin cement with a universal bonding agent.
-  **Easy Cleanup** - Cleanup is a cinch with One-Peel™ removal and with Maxcem Elite Chroma's Cleanup Indicator which dispenses pink before fading at the gel state, indicating the optimal time to clean up excess cement.
-  **Universal Compatibility** - Simplify with one cement for your indirect restoration needs. Maxcem Elite and Maxcem Elite Chroma Universal Resin Cement are compatible with all* dental substrates and indirect restorations including anterior or posterior ceramics, PFMs and metal restorations, and CAD/CAM materials.



SimpliShade™ Bulk Fill

One Shade, one Step

-  **One-Shade System** - Using one shade that blends seamlessly with all 16 VITA® classical shades, SimpliShade Bulk Fill saves you time selecting shades and managing inventory.
-  **Chameleon Effect** - Special nanofiller architecture mimics the natural enamel structure and interacts with light for virtually indistinguishable blending.
-  **One Step** - Eliminate the need for a capping layer, complete restorations with a SingleFill™ composite. With its enhanced curing mechanism, SimpliShade Bulk Fill makes it easy to achieve polymerization depths up to 5mm in a single step.



Temp-Bond™

Temporary attachment. Permanent reliability.

-  **High bond strength** - Designed to prevent leakage.
-  **Easy to remove** - Minimises potential damage to the prep or temporary.
-  **Flexible delivery choices** - Temp-Bond Original, Temp-Bond Clear and Temp-Bond NE are available in automix syringe. Moreover, Temp-Bond and Temp-Bond NE are also available in Unidose® and tubes.
-  **Flows and mixes easily** - Delivers optimal consistency for solid, complete seating of restorations.



SimpliShade™ Bulk Fill Flow

Simplicity in a Single-Shade Bulk Fill Flowable

-  **One-Shade System** - SimpliShade Bulk Fill Flow offers single-shade simplicity with the durability and depth of cure you demand.
-  **Save Time** - Simplify your workflow with just one step to cure up to 5 mm, before adding a capping layer.
-  **Go with the Flow** - Fill every area of the tooth's anatomy with a self-leveling bulk fill flowable that easily adapts to the cavity floor.



Demi™ Plus

LED Light Curing System



True 5-Second Cure - Consistent output of Periodic Level Shifting (PLS) delivers a true 5-second cure for composite shades A3 and lighter.



Lightweight, Ergonomic Design - Multiple holding options can reduce hand and arm stress.



OptiDisc™

High performing polishing and contouring disc.



Flexible and precise - The shorter madrel and the disc's bendability allow you to reach posterior and other restricted areas with ease. The transparent material gives you great visibility while working.



Safe - OptiDisc's active side is easy to recognise and the patented mandrel locks below the surface of the disc, protecting your restoration against scarring and scratching.



4 colours, 4 stages of abrasiveness - The OptiDisc system offers the choice of four grit level discs for contouring, finishing and high gloss polishing. It is suitable for composite, glass ionomers, amalgams, semiprecious and precious metals. Discs are also available in three different diameter sizes depending on your needs: 9.6mm, 12.6mm and 15.9mm.



Identoflex™ Diamond Ceramic Polishers

Chairside diamond polishers for ceramics.



Safe - Thanks to the high polishing power of diamond, lower revs can be used. Overheating of teeth, which can damage the pulp and filling materials, is thus avoided.



Saving resources - Twin-zone technology takes account of typical wear patterns and reduces wear of diamond grit. The active zone is easily identifiable due to colour-coding.



Occlubrush™

One-step polishing brush system.



Designed for all surfaces - Non-destructive polishing of occlusal, concave and difficult to reach surfaces



Easy and efficient - Polishing particles already embedded in the bristles, so no paste is necessary.



OptiShine™

Concave-shaped polishing brush.



Universal use - For restorations made of composite, compomer, resin-modified glass ionomer or ceramics, as well as for prophylactic treatment.



Efficient in practice - The polishing effect is created by polishing particles embedded in the bristles (silicon carbide), so no paste is necessary.



Multiple use - Autoclavable at 134°C for at least 3 minutes. No effect on the polishing performance.



Opti1Step™

Composite polishing system.



Universal use - Different shapes available to access different areas and angles. Allows variable abrasion by changing pressure.



Easy and fast technique - No paste needed, for all types of composite restorations, designed for either wet or dry applications.



Multiple use - Autoclavable at 134°C for at least 3 minutes.

ORDERING INFORMATION

Demi™ Plus Curing Light

910860-2 Demi Plus LED Light Curing System (220 V)
Contains: Handpiece, 8 mm Turbo Extended Light Guide,
Lithium Ion Battery Pack, Battery Charger, Power Supply,
Protective Light Shield.

Identoflex™ Diamond Ceramic Polishers

ID 5721/6 6 x Flame

ID 5781/6 6 x Cup

Maxcem Elite™ Chroma Universal Resin Cement

36371 Standard Kit

36372 Clear Refill (2 dual syringes)

36373 White Refill (2 dual syringes)

36374 Yellow Refill (2 dual syringes)

Maxcem Elite™ Universal Resin Cement

34055 Standard Kit

34056 Value Kit

34058 Clear Refill (2 dual syringes)

34059 White Refill (2 dual syringes)

34060 White Opaque Refill (2 dual syringes)

34061 Yellow Refill (2 dual syringes)

Occlubrush™

2510 10x Cup

2511 10x small Cup

2512 10x Point

2520 Set with 5 of each

Opti1Step™ Polisher

8001 12 x Flame

8002 12 x Cup

8003 12 x Disc

OptiBond™ Universal

36517 Bottle Kit

36518 Unidose Refills (pkg of 100)

36519 Bottle Refill (5ml)

26909 Gel Etchant Dispensing Tips (pkg of 30)

31297 Gel Etchant Refill Pack

OptiDisc™

4200 General Assorted Kit

4190 Assorted Kit

4188 Mini Kit

Refills 100 pcs (except Extra-Course, 80 pcs)

4197 9.6mm – Extra-Coarse **4185** 12.6mm – Fine

4181 9.6mm – Coarse/Medium **4186** 12.6mm – Extra-Fine

4182 9.6mm – Fine **4199** 15.9mm – Extra-Coarse

4183 9.6mm – Extra-Fine **4191** 15.9mm – Coarse/Medium

4198 12.6mm – Extra-Coarse **4192** 15.9mm – Fine

4184 12.6mm – Coarse/Medium **4193** 15.9mm – Extra-Fine

OptiShine™

2513 3 Brushes

2514 10 Brushes

SimpliShade™ Bulk Fill

37211 Syringe

37212 Unidose 20 Pack

SimpliShade™ Bulk Fill Flow

37217 Syringe

37218 Unidose 20 Pack

SimpliShade™ Universal Composite

37001 Syringe, Light

37002 Syringe, Medium

37003 Syringe, Dark

37004 Syringe, Bleach White

37005 Syringe, Universal Opaque

37006 Unidose 20 Pack, Light

37007 Unidose 20 Pack, Medium

37008 Unidose 10 Pack, Dark

37009 Unidose 10 Pack, Bleach White

37010 Unidose 10 Pack, Universal Opaque

Take 1™ Advanced™

34963 Fast Set

34962 Regular Set

Temp-Bond™

60058 Introductory Pack (with Modifier)

61087 Standard Pack (without Modifier)

33684 2 x Automix Syringe

31377 Unidose 50 x 2.4g

Temp-Bond™ Clear

33351EU 2 x Automix Syringe

33361 50 x Automix Tip Refill

Temp-Bond™ NE

60256 Standard Pack

33685 2 x Automix Syringe

31973 Unidose 50 x 2.4g

TempSpan™ Temporary Crown & Bridge Material

N69AA Shade A0

N69AE Shade A3.5

N69AB Shade A1

N69AF Shade B1

N69AC Shade A2

N69AG Shade C2

N69AD Shade A3

Contact your Kerr Rep for a complete list of rotary SKUs



For More Information Scan the QR Code or
Visit to.kerrdental.com/en/simplicity