



- ZFO has an MFD of 1 mm
- ZFO Triangular cross section has 3 positive cutting edges vs 2 for WOG
- ZFO chip space area is considerably greater than WOG, reducing apical debris extrusion significantly
- 'R' Phase proprietary Variable HT based on files Mass volume

kerrdental.com

©2026 Kerr Corporation. All trademarks are property of Kerr Corporation. All Rights Reserved. MKT-26-0505 Rev1



ZenFlex™ ONE Reciprocating File System



Size	ISO size 20 ● Small	ISO size 25 ● Primary	ISO size 35 ● Medium	ISO size 45 ● Large
Taper at the tip	6%	6%	6%	6%
Lengths	21 mm / 25 mm / 31 mm			

ZenFlex™ ONE Guttha Percha Points



Size	ISO size 20 ●	ISO size 25 ●	ISO size 35 ●	ISO size 45 ●
Taper at the tip	6%	6%	6%	6%
Lengths	30mm			

ZenFlex™ ONE Paper Points



Size	One Size
	Fits to all (Small, Primary, Medium, Large)
Lengths	30mm

MKT-26-0505_Rev1